



# DIAMOND

## DIAMOND TIP DERMABRASION



***A medical & cosmetic alternative for skin treatments.***

Diamond tip Dermabrasion is a system by which the skin is slightly suctioned by negative pressure thus allowing to be abraded at different levels by circular heads of varied micronage.

**SKIN AGING • WRINKLES AND FURROWS  
BLACKHEADS • HYPERCHROMIA  
STRETCH MARKS • SCARS**



**abrasion  
tips**

### SPECIFICATIONS

Design and size intended for market requirements  
Weight: 4 K  
Dimensions: Length 280 mm, width 230 mm, height 125 mm  
Vacuum Output: -0.8 bar  
Vacuum Gauge - Background Scale: -1,0 bar  
Power: 90 W

### EASY TO CARRY

**DIAMOND TIPS:**  
**22 mm Diameter:** 180 UM and 100 UM  
**19 mm Diameter:** 220 UM and 180 UM  
**15 mm Diameter:** 250 UM, 200UM and 100UM  
**8 mm Diameter:** 200 UM  
**6 mm Diameter:** 250 UM

## EQUIPMENT

**Diamond** features nine different diamond tip options according to the desired depth of decortication.


Aged skin, dead cells and other debris are removed by abrasion, in a way that superficial tissues are renewed and maintained in optimal conditions. The use of diamond tip dermabrasion allows therapists to carry out their treatments in a controlled and precise way, free from residual material.

**It is an essential tool for professionals who work on the skin, as it has 9 diamond tips of different micronage and diameter.**

The vacuum gauge quantifies the pressure applied on the patient during the procedure and this helps keeping track of the treatment.

### Procedure

**Diamond** can be used in various facial and body treatments. The finer abrasive tip achieves a light polishing of the skin. The medium are ideal for more aggressive facial exfoliation and the coarse tips work well when exfoliating chest, hands, arms, elbows and other body regions. The vacuum level can be adjusted to simultaneously increase blood flow and achieve oxygenation of the skin.



## DIAMOND TIPS



**22 mm Diameter:** Body, especially for acne, knee, elbow and areas with stretch marks.



**19 mm Diameter:** Face, especially for large stains, pores, premature aging due to sun exposure, etc.



**15 mm Diameter:** Scars, fine lines and various imperfections. Soft Dermabrasion.



**8mm y 6mm Diameter:** Eyes and both sides of the nose. Wrinkles, pimples and blackheads, acne.